

## MEDICAL RAIL CLAMPS FOR EQUIPMENT RAILS

EU DIN standard 10x25 mm



**Rail clamp, narrow model, with one ball lock and T-slot**

Item no.: **JD1460000**  
 Material: **Aluminium (Natural Anodized)**  
 Dimensions: **Height x Width x Depth: 61,5 mm x 64 mm x 28,6 mm**  
 Weight: **132 grams**  
 MOQ: **5 pcs.**



**Rail clamp narrow model, with one ball lock**

Item no.: **JD1470000**  
 Material: **Aluminium (Natural Anodized)**  
 Dimensions: **Height x Width x Depth: 61,5 mm x 64 mm x 28,6 mm**  
 Weight: **168 grams**  
 MOQ: **5 pcs.**



**Rail clamp, narrow model, with one ball lock & three countersunk Ø6,6 holes**

Item no.: **JD1470300**  
 Material: **Aluminium (Natural Anodized)**  
 Dimensions: **Height x Width x Depth: 61,5 mm x 64 mm x 28,6 mm**  
 Weight: **161 grams**  
 MOQ: **5 pcs.**

“Prepared for suspension.”



**Rail clamp, wide model, locked using two hex allen screws**

Item no.: **JD1430000**  
 Material: **Aluminium (Natural Anodized)**  
 Dimensions: **Height x Width x Depth: 61,5mm x 90 mm x 28,6 mm**  
 Weight: **308 grams**  
 MOQ: **5 pcs.**

## RAIL CLAMPS

EU DIN standard for 10x25 mm clamp rails



**Rail clamp, wide model, locked using two hex allen screws. It has three countersunk holes Ø6,6 mm.**

“Prepared for suspension.”

**Rail clamp, wide model, locked with two ball locks.**

**Rail clamp, wide model, locked with two ball locks. It has three countersunk holes Ø6,6 mm.**

“Prepared for suspension.”

**Rail clamp, wide model, locked with one ball lock.**

Item no.: **JD1430300**  
 Material: **Aluminium (Natural Anodized)**  
 Dimensions: **Height x Width x Depth: 61,5 mm x 90 mm x 28,6 mm**  
 Weight: **303 grams**  
 MOQ: **5 pcs.**


Item no.: **JD1440000**  
 Material: **Aluminium (Natural Anodized)**  
 Dimensions: **Height x Width x Depth: 61,5 mm x 110 mm x 28,6 mm**  
 Weight: **328 grams**  
 MOQ: **5 pcs.**

Item no.: **JD1440300**  
 Material: **Aluminium (Natural Anodized)**  
 Dimensions: **Height x Width x Depth: 61,5 mm x 110 mm x 28,6 mm**  
 Weight: **318 grams**  
 MOQ: **5 pcs.**

Item no.: **JD1450000**  
 Material: **Aluminium (Natural Anodized)**  
 Dimensions: **Height x Width x Depth: 61,5mm x 110 mm x 28,6 mm**  
 Weight: **325 grams**  
 MOQ: **5 pcs.**

## RAIL CLAMPS

EU DIN standard for 10x25 mm clamp rails



**Rail clamp, wide model, locked with one ball lock. It has three countersunk holes Ø6,6 mm.**

“Prepared for suspension.”

**Rail clamp, narrow model, locked with two socket screws. It has three countersunk Ø6,6 mm holes.**

“Prepared for suspension.”

**Rail clamp, wide model, locked using two hex allen screws & adapter fitting Ø18 mm**

Three bush sizes are available.  
 Handles can be purchased.

**Rail clamp, wide model, 1 pcs. ball lock & adapter fitting Ø18 mm**

Three bush sizes are available.  
 Handles can be purchased.

Item no.: **JD1450300**  
 Material: **Aluminium (Natural Anodized)**  
 Dimensions: **Height x Width x Depth: 61,5 mm x 110 mm x 28,6 mm**  
 Weight: **317 grams**  
 MOQ: **5 pcs.**

Item no.: **JD0450300**  
 Material: **Aluminium (Natural Anodized)**  
 Dimensions: **Height x Width x Depth: 61,5 mm x 45 mm x 28,6 mm**  
 Weight: **151 grams**  
 MOQ: **5 pcs.**

Item no.: **JD1430318**  
 Material: **Aluminium (Natural Anodized)**  
 Dimensions: **Lenght x Width x Depth: 61,5 mm x 90 mm x 64,2 mm**  
 Weight: **468 grams**  
 MOQ: **2 pcs.**

Item no.: **JD1450318**  
 Material: **Aluminium (Natural Anodized)**  
 Dimensions: **Lenght x Width x Depth: 61,5 mm x 110 mm x 64,2 mm**  
 Weight: **501 grams**  
 MOQ: **2 pcs.**

## LAMP RAIL CLAMPS - EXTRA SOLID - FOR LONG AND HEAVY LAMPS

EU DIN STANDARD FOR 10X25



## SUSPENSION RAIL CLAMPS FOR DROP STANDS, IT & MEDICAL EQUIPMENT

EU DIN STANDARD FOR 10X25 MM

**Rail clamp, wide model, 2 pcs. ball locks & adapter fitting Ø18 mm**

**Rail clamp, narrow model, 1 pcs. ball locks & adapter fitting Ø18 mm**

**Rail clamp, wide model, locked using two hex allen screws & adapter fitting Ø18 mm**

**Rail clamp, wide model, with a one-ball lock, adapter fitting & brass bush Ø20 mm. hole**

Three bush sizes are available.  
Handles can be purchased

Three bush sizes available .  
Handles can be purchased

Handles can be purchased

Handles can be purchased

Item no.: **JD1440318**  
 Material: **Aluminium (Natural Anodized)**  
 Dimensions: **Length x Width x Depth: 61,5 mm x 110 mm x 64,2 mm**  
 Weight: **502 grams**  
 MOQ: **2 pcs.**

Item no.: **JD1470318**  
 Material: **Aluminium (Natural Anodized)**  
 Dimensions: **Length x Width x Depth: 61,5 mm x 64 mm x 64,2 mm**  
 Weight: **343 grams**  
 MOQ: **4 pcs.**

Item no.: **JD1430320**  
 Material: **Aluminium (Natural Anodized)**  
 Dimensions: **Length x Width x Depth: 61,5 mm x 90 mm x 60,6 mm**  
 Weight: **536 grams**  
 MOQ: **2 pcs.**

Item no.: **JD1450320**  
 Material: **Aluminium (Natural Anodized)**  
 Dimensions: **Length x Width x Depth: 61,5mm x 110 mm x 64,2 mm**  
 Weight: **550 grams**  
 MOQ: **2 pcs.**

## ACCESSORIES FOR VARIOUS LAMP RAIL CLAMPS

### POM bushing with bushing hole $\varnothing$ 9.6mm

Item no.: **JD128009**  
Material: **Plastic (White)**  
Weight: **12 grams**  
MOQ: **10 pcs.**



### Brass bushing with bushing hole $\varnothing$ 20 mm

Item no.: **JD220005**  
Material: **Brass**  
Weight: **88 grams**  
MOQ: **10 pcs.**



### POM bushing with bushing hole $\varnothing$ 13mm

Item no.: **JD128001**  
Material: **Plastic (White)**  
Weight: **8 grams**  
MOQ: **10 pcs.**



### Handle with M6 x 20mm thread

Item no.: **JD940000**  
Material: **Stainless Steel (AISI 304)**  
Weight: **44 grams**  
MOQ: **5 pcs.**



### POM bushing with bushing hole $\varnothing$ 16mm

Item no.: **JD1270001**  
Material: **Plastic (White)**  
Weight: **3,6 grams**  
MOQ: **10 pcs.**



### Handle with M10 x 40mm thread

Item no.: **JD460000**  
Material: **Stainless Steel (AISI 304)**  
Weight: **120 grams**  
MOQ: **2 pcs.**



## MULTIBRACKET - ADJUSTABLE - 16-41MM



**Multibracket - The fitting can be connected anywhere. Fits Ø16-41 mm.**

This multibracket makes it possible to assemble various equipment together.

Item no.: **JD1580000**  
 Material: **Aluminium (Natural Anodized), SS 304, POM og Nylon (PGA)**  
 Dimensions: **Lenght x Width x Depth: 36 mm x 51 mm x 134/96 mm**  
 Weight: **209 grams**  
 MOQ: **4 pcs.**



**Multibracket with 2xØ6 clearance holes. Fits up to. Ø16-41 mm.**

This multibracket makes it possible to assemble various equipment together.

Item no.: **JD1580005**  
 Material: **Aluminium (Natural Anodized), SS 304, POM og Nylon (PGA)**  
 Dimensions: **Lenght x Width x Depth: 36 mm x 51 mm x 134/96 mm**  
 Weight: **204 grams**  
 MOQ: **4 pcs.**



**Double multibracket for mounting of for example bed & drop stand.**

Includes 2 pcs. Unbrako bolt screws.

Item no.: **JD1580605**  
 Material: **Aluminium (Natural Anodized), SS 304, POM og Nylon (PGA)**  
 Weight: **421 grams**  
 MOQ: **2 pcs.**



**Multibracket with T-slot fitting fits from Ø16-41 mm. Fitting can be placed on table, trolley, drop stand, DIN rail and columns.**

Includes 2 pcs. hex allen screws.

Item no.: **JD15800163**  
 Material: **Aluminium (Natural Anodized), SS 304, POM og Nylon (PGA)**  
 Dimensions: **Lenght x Width x Depth: 61,5mm x 64 mm x 60,5 mm**  
 Weight: **245 grams**  
 MOQ: **4 pcs.**

## MULTIBRACKET - ADJUSTABLE - 16-41MM - CONTINUED



**Multibracket with a fixing device, Ø18mm hole.  
Fits Ø16-41 mm.**

Three bush sizes available .  
Handle can be purchased separately.

Item no.: **JD1580318**  
 Material: **Aluminium (Natural Anodized),  
SS 304, POM, Nylon (PA6)**  
 Dimensions: **Lenght x Width x Depth:  
36 mm x 51 mm x 134/96 mm**  
 Weight: **288 grams**  
 MOQ: **2 pcs.**

**ACCESSORIES:** Handle can be purchased.



**Multibracket with a fixing device and brass bushing.  
Fits up to. Ø16-41 mm.**

Handle can be purchased separately.

Item no.: **JD1580320**  
 Material: **Aluminium (Natural Anodized),  
SS 304, Brass, POM, Nylon (PA6)**  
 Dimensions: **Lenght x Width x Depth:  
36 mm x 51 mm x 134/96 mm**  
 Weight: **440 grams**  
 MOQ: **2 pcs.**

**ACCESSORIES:** Handle and T-slot fitting can be purchased.



**T-slot fitting for multi-bracket**

The holes are only suitable for multibracket.

Item no.: **JD1630001**  
 Material: **Aluminium (Natural Anodized)**  
 Dimensions: **Height x Width x Depth:  
60 mm x 36,2 mm x 12,5 mm**  
 Weight: **36 grams**  
 MOQ: **1 Bag med 20 pcs.**

## FITTINGS SUITABLE FOR GCX RAILS



**GCX rail fitting, prepared with mounting holes**

Item no.: **JD82062506**  
 Material: **Aluminium (Natural Anodized)**  
 Colour: **Natural**  
 Dimensions: **Height x Width x Depth: 82 mm x 63 mm x 6 mm**  
**83 grams**  
 Weight:  
 MOQ: **4 pcs.**



**GCX rail fitting with a T-slot.**

Item no.: **JD820625163**  
 Material: **Aluminium (Natural Anodized)**  
 Colour: **Natural**  
 Dimensions: **Height x Width x Depth: 82 mm x 36 mm x 18,5 mm**  
**123 grams**  
 Weight:  
 MOQ: **4 pcs.**



**GCX rail fitting with suspension housing. Ø20 mm bushing hole.**

Item no.: **JD82062520**  
 Material: **Aluminium (Natural Anodized) & Brass**  
 Colour: **Natural**  
 Dimensions: **Height x Width x Depth: 82 mm x 63 mm x 42 mm**  
**313 grams**  
 Weight:  
 MOQ: **2 pcs.**

**ACCESSORIES:** Handle can be purchased.



*POM bushing not included as standard*

**GCX rail fitting with suspension housing Ø18 mm hole.**

Optional plastic bushing Ø9,6/13 or 16 mm.

Item no.: **JD82062518**  
 Material: **Aluminium (Natural Anodized) & Plastic**  
 Colour: **Natural**  
 Dimensions: **Height x Width x Depth: 82 mm x 63 mm x 42 mm**  
**270 grams**  
 Weight:  
 MOQ: **2 pcs.**

**ACCESSORIES:** Handle can be purchased.



## FITTINGS SUITABLE FOR GCX RAILS



For flat screen / monitor holder



For catheter holder



For light fixture



For sharpsafe holder

## FITTINGS SUITABLE & GCX RAILS



For monitor holder



For wire basket




Fitting with suspension housing



JD820625163  
For catheter holder

## GLOVE, WET WIPE & SANITARY FACEMASK HOLDER



With “New HandsUp”, staff can finally see the glove sizes, which are always written on the sides of the boxes.

Fits 99 percent of people glove boxes.

Fits Wet Wipes (the large MAXI)

Fits sanitary towels and mask boxes

“New HandsUp” can be mounted in both horizontal and vertical position.

Mounts to wall fitting, T-slot or wall rail

Autoclavable

Item no.:  
Material:

Surface:  
Dimensions:

Weight:  
MOQ:

**JD8000000**  
**Stainless Steel (AISI 304) -**  
**4mm spring steel**  
**Electropolished**  
**Depth x Height x Width**  
**105 mm x 148 mm x 260 mm**  
**362 grams**  
**4 pcs.**



## GLOVE BOX HOLDER



**Glove holder for 3 glove boxes wall mounted**

Suitable for glove boxes, face wipes, aprons and handkerchiefs, etc.

Item no.: **JD1100000**  
 Material: **Stainless Steel (AISI 304)**  
 Dimensions: **Depth x Height x Width (outside)**  
**93 mm x 406 mm x 285 mm**  
 Opening: **183 mm**  
 Weight: **1.854 grams**  
 MOQ: **4 pcs.**



**Glove holder for 4 glove boxes wall mounted**

Suitable for glove boxes, face wipes, aprons and handkerchiefs, etc.

Item no.: **JD1100004**  
 Material: **Stainless Steel (AISI 304)**  
 Dimensions: **Depth x Height x Width (outside)**  
**93 mm x 520 mm x 285 mm**  
 Opening: **183 mm**  
 Weight: **2.346 grams**  
 MOQ: **4 pcs.**

# SHARPS CONTAINER HOLDER - ROUND FOR E.G. "Grathwol"



## Sharps Container 0,7 & 1,5 L

Sharps container 0,7 & 1,5 liter  
AP Medical DISPO Line

Yellow or red lid

### Accessories:

- Wall bracket with T-slot
- Rail clamp with T-slot (DK)
- Rail clamp with T-slot (EU)
- Multibracket with T-slot
- GCX rail bracket with T-slot



Sharps container holder til e.g.  
AP Medical 0,7 & 1,5 liter Sharps Containers

|                  |                                   |
|------------------|-----------------------------------|
| Item no.:        | <b>JD2550000</b>                  |
| Inner Diameter:  | <b>Ca. Ø100 mm</b>                |
| Material:        | <b>Stainless Steel (AISI 304)</b> |
| Surface:         | <b>Electropolished</b>            |
| Height:          | <b>50 mm</b>                      |
| Depth:           | <b>143 mm (outer dimension)</b>   |
| Thread Diameter: | <b>8 mm</b>                       |
| Weight:          | <b>189 grams</b>                  |
| MOQ:             | <b>5 pcs.</b>                     |

## SHARPS CONTAINER HOLDER - ROUND FOR E.G. "AP MEDICAL"



### Sharps Container 2 Liter

Sharps container 2 litres  
AP Medical DISPO Line

Yellow or red lid



Sharps container holder for  
1.0-1.3-1.5-2.0-3.0 litre Sharps Containers

Fits perfectly to Klinion EasyCare "Mediq /UD-SON" 1,3 & 2L and 2 liter AP Medical DISPO Line boxes. Yellow or red lid.

|                  |                                   |
|------------------|-----------------------------------|
| Item no.:        | <b>JD1490000</b>                  |
| Inner Diameter:  | <b>Ø 134 mm</b>                   |
| Material:        | <b>Stainless Steel (AISI 304)</b> |
| Surface:         | <b>Electropolished</b>            |
| Height:          | <b>50 mm</b>                      |
| Depth:           | <b>172 mm (outer dimension)</b>   |
| Thread Diameter: | <b>8 mm</b>                       |
| Weight:          | <b>220 grams</b>                  |
| MOQ:             | <b>5 pcs.</b>                     |

## SHARPS CONTAINER HOLDER - OVAL FOR E.G. "AP MEDICAL"



### Sharps Container 2, 4 & 7 Liter

Sharps container 2 litres  
AP Medical DISPO Line

Yellow or red lid

**Accessories:**  
 Wall bracket with T-slot  
 Rail clamp with T-slot (DK)  
 Rail clamp with T-slot (EU)  
 Multibracket with T-slot  
 GCX rail bracket with T-slot

Sharps container holder for e.g. AP Medical 2, 4  
& 7 litre oval boxes

Fits perfectly to Klinion EasyCare "Mediq /UD-SON" AP Medical DISPO Line boxes.  
Yellow or red lid.

|                  |                                   |
|------------------|-----------------------------------|
| Item no.:        | <b>JD1590000</b>                  |
| Inner Diameter:  | <b>133 mm x 133 mm</b>            |
| Material:        | <b>Stainless Steel (AISI 304)</b> |
| Surface:         | <b>Electropolished</b>            |
| Height:          | <b>50 mm</b>                      |
| Depth:           | <b>175 mm (outer dimension)</b>   |
| Thread Diameter: | <b>8 mm</b>                       |
| Weight:          | <b>248 grams</b>                  |
| MOQ:             | <b>5 pcs.</b>                     |

## SHARPS CONTAINER HOLDER - OVAL

For E.G. "AP Medical"



### Sharps Container 3 Liter

Sharps container 2 liter AP Medical DISPO Line

Yellow or red lid

#### Accessories:

Wall bracket with T-slot  
 Rail clamp with T-slot (DK)  
 Rail clamp with T-slot (EU)  
 Multibracket with T-slot  
 GCX rail bracket with T-slot

Sharps container holder for e.g. AP Medical  
 3 liter, oval container

Item no.: **JD1690000**  
 Inner Diameter: **143 mm x 143 mm**  
 Material: **Stainless Steel (AISI 304)**  
 Surface: **Electropolished**  
 Height: **50 mm**  
 Depth: **175 mm (outer dimension)**  
 Thread Diameter: **8 mm**  
 Weight: **319 grams**  
 MOQ: **5 pcs.**

## SHARPS CONTAINER HOLDER - SQUARED

For E.G. "Sharpsafe"



### Sharps Container 4 & 7 Liter

Sharps containers 4 / 7 liter  
 "Sharpsafe" container

Yellow or red lid

Sharps container holder for e.g. "Sharpsafe"  
 4 & 7 liter squared container

Item no.: **JD1540000**  
 Inner Diameter: **156 mm x 156 mm**  
 Material: **Stainless Steel (AISI 304)**  
 Surface: **Electropolished**  
 Height: **50 mm**  
 Depth: **208 mm (outer dimension)**  
 Thread Diameter: **8 mm**  
 Weight: **280 grams**  
 MOQ: **5 pcs.**



## SHARPS CONTAINER HOLDER - SQUARED

For E.G. "Sharpsafe"



### Sharps Container 7,5 Liter

Sharps container "Sharpsafe"  
7.5 liter QUIVER

Yellow or red lid

Sharps container holder for e.g.  
"Sharpsafe Quiver"  
7,5 liter, squared container (long model)

Item no.: **JD2790000**  
 Inner Diameter: **101 mm x 147 mm**  
 Material: **Stainless Steel (AISI 304)**  
 Surface: **Electropolished**  
 Height: **50 mm**  
 Depth: **153 mm (outer dimension)**  
 Thread Diameter: **8 mm**  
 Weight: **247 grams**  
 MOQ: **5 pcs.**

**Accessories:**  
 Wall bracket with T-slot  
 Rail clamp with T-slot (DK)  
 Rail clamp with T-slot (EU)  
 Multibracket with T-slot  
 GCX rail bracket with T-spur

## SHARPS CONTAINER HOLDER - SQUARED

For E.G. "Europlásticos"



### Sharps Container 1 Liter

Sharps container for 1 liter boks  
from Europlásticos

Yellow or red lid

Sharps container holder for e.g. "Sharpsafe"  
1 liter squared container

Item no.: **JD1700000**  
 Inner Diameter: **96 mm x 96 mm**  
 Material: **Stainless Steel (AISI 304)**  
 Surface: **Electropolished**  
 Height: **50 mm**  
 Depth: **208 mm (outer dimension)**  
 Thread Diameter: **8 mm**  
 Weight: **200 grams**  
 MOQ: **5 pcs.**





## SHARPS CONTAINER HOLDER - SQUARED

Other Assortment - e.g. "Safebox"



### Sharps Container 4L, 6L & 11 Liter

Sharps Container "Safebox" 4L, 6L & 11L

Yellow or red lid

#### Accessories:

Wall bracket with T-slot  
 Rail clamp with T-slot (DK)  
 Rail clamp with T-slot (EU)  
 Multibracket with T-slot  
 GCX rail bracket with T-slot



Sharps Container Holder for e.g.  
 "Sharpsafe Quiver"  
 7,5 liter, squared container (long model)

Item no.: **JD2790000**  
 Inner Diameter: **101 mm x 147 mm**  
 Material: **Stainless Steel (AISI 304)**  
 Surface: **Electropolished**  
 Height: **50 mm**  
 Depth: **153 mm (outer dimension)**  
 Thread Diameter: **8 mm**  
 Weight: **247 grams**  
 MOQ: **5 pcs.**

## SHARPS CONTAINER HOLDER - ROUND

Other Assortment



Sharps Container Holder for  
 0,5- & 1 liter round container

Fits e.g. Dispo 0,5 litre box from Grathwol.

Item no.: **JD1500000**  
 Inner Diameter: **Ø73 mm**  
 Material: **Stainless Steel (AISI 304)**  
 Surface: **Electropolished**  
 Height: **50 mm**  
 Depth: **112 mm (outer dimension)**  
 Thread Diameter: **8 mm**  
 Weight: **150 grams**  
 MOQ: **5 pcs.**

Sharps Container Holder for  
 0.6 litre - round container

Fits e.g. Multisafe 1 l. & USON 500 ml.

Item no.: **JD1510000**  
 Inner Diameter: **Ø89 mm**  
 Material: **Stainless Steel (AISI 304)**  
 Surface: **Electropolished**  
 Height: **50 mm**  
 Depth: **134 mm (outer dimension)**  
 Thread Diameter: **8 mm**  
 Weight: **171 grams**  
 MOQ: **5 pcs.**

## SHARPS CONTAINER HOLDER - ROUND

Other Assortment



**Sharps Container Holder for e.g. Mediq /USON  
5 & 7 litre round container**

Item no.: **JD2640000**  
 Inner Diameter: **Ø200 mm**  
 Material: **Stainless Steel (AISI 304)**  
 Surface: **Electropolished**  
 Height: **50 mm**  
 Depth: **244 mm (outer dimension)**  
 Thread Diameter: **8 mm**  
 Weight: **286 grams**  
 MOQ: **5 pcs.**

**Sharps Container Holder for e.g. USON 11 litres  
Jokey container from Mediq**

Item no.: **JD2660000**  
 Inner Diameter: **Ø251 mm**  
 Material: **Stainless Steel (AISI 304)**  
 Surface: **Electropolished**  
 Height: **50 mm**  
 Depth: **295 mm (outer dimension)**  
 Thread Diameter: **8 mm**  
 Weight: **345 grams**  
 MOQ: **5 pcs.**


**Sharps Container Holder for e.g. USON 21 litre  
and Jokey container from Mediq**

Item no.: **JD2670000**  
 Inner Diameter: **Ø280 mm**  
 Material: **Stainless Steel (AISI 304)**  
 Surface: **Electropolished**  
 Height: **50 mm**  
 Depth: **324 mm (outer dimension)**  
 Thread Diameter: **8 mm**  
 Weight: **368 grams**  
 MOQ: **5 pcs.**

**Accessories:**  
 Wall bracket with T-slot  
 Rail clamp with T-slot (DK)  
 Rail clamp with T-slot (EU)  
 Multibracket with T-slot  
 GCX rail bracket with T-slot

## WIRE BASKETS & WIRE PRODUCTS

Stainless steel



**Wire basket 6 litres, conical  
with hook for T-slot**

**Wire basket 7 litres, conical  
with hook for T-slot**

**Wire basket 15 litres, conical  
with hook for T-slot**

Item no.: **JD1610000**  
 Dimensions: **Width x Depth x Height  
231mm x 161mm x 242mm  
Stainless Steel (AISI 304)**  
 Material: **Stainless Steel (AISI 304)**  
 Surface: **Electropolished**  
 Weight: **853 grams**  
 Volumen: **6 Liter**  
 MOQ: **5 pcs.**

**Accessories for all 3 types of wire baskets:**  
 Wall bracket with T-slot  
 Rail clamp with T-slot (DK)  
 Rail clamp with T-slot (EU)  
 Multibracket with T-slot  
 GCX rail bracket with T-spur

Item no.: **JD1610101**  
 Dimensions: **Width x Depth x Height  
623mm x 217mm x 132mm  
Stainless Steel (AISI 304)**  
 Material: **Stainless Steel (AISI 304)**  
 Surface: **Electropolished**  
 Weight: **667 grams**  
 Volumen: **7 Liter**  
 MOQ: **5 pcs.**

Item no.: **JD1610201**  
 Dimensions: **Width x Depth x Height  
337mm x 233mm x 267mm  
Stainless Steel (AISI 304)**  
 Material: **Stainless Steel (AISI 304)**  
 Surface: **Electropolished**  
 Weight: **959 grams**  
 Volumen: **15 Liter**  
 MOQ: **5 pcs.**

## WIRE BASKETS & WIRE ACCESSORIES



**Wire basket 25 litres, conical**

Item no.: **JD1610300**  
 Dimensions: **Width x Depth x Height**  
**348mm x 248mm x 400mm**  
**Stainless Steel (AISI 304)**  
 Surface: **Electropolished**  
 Weight: **1.241 grams**  
 Volumen: **25 Liter**  
 MOQ: **5 pcs.**



**Wall bracket for hanging 25 litre wire basket**

Item no.: **JD2810000**  
 Dimensions: **Width x Depth x Height**  
**23,5mm (above) 13,8mm**  
**(below) x 21mm x 53,5mm**  
**Stainless Steel (AISI 304)**  
 Material: **Electropolished**  
 Surface: **Electropolished**  
 Weight: **27 grams**  
 Hole size: **Ø4,25 mm**  
 MOQ: **40 pcs. / bag**

## WALL BRACKETS FOR BED RACKS

Max Ø38 mm tube



**Mask & tubing hook, with a guide for T-slot**

Item no.: **JD1620000**  
 Dimensions: **Width x Depth x Height**  
**30mm x 76mm x 100mm**  
**Stainless Steel (AISI 304)**  
 Surface: **Electropolished**  
 Weight: **80 grams**  
 MOQ: **10 pcs.**



**Wall bracket for hanging bed equipment**

Enables safe suspension of various hospital & care items.

Item no.: **JD1464538**  
 Dimensions: **Width (above) or Width (below) x Depth x Height**  
**40mm (above) 20mm (below) x 50mm x 93 mm**  
**Stainless Steel (AISI 304)**  
 Material: **Electropolished**  
 Surface: **Electropolished**  
 Weight: **117 grams**  
 Hole size: **Ø6 mm**  
 MOQ: **10 pcs.**

## BRACKET FOR SUCTION CANISTER & CATHETER CONTAINER



**Bracket Ø100, for 2 Litre  
suction canister, Ø100 mm. hole**



**Catheter container  
with guide for T-slot**

Item no.: **JD1000000**  
Dimensions: **Ø100 mm / Ø 8mm Thread**  
Material: **Stainless Steel (AISI 304)**  
Surface: **Electropolished**  
Weight: **189 grams**  
MOQ: **5 pcs.**

Item no.: **JD2390000**  
Height: **400 mm**  
Material: **Stainless Steel (AISI 304)**  
Surface: **Electropolished**  
Weight: **1.000 grams**  
Volumen: **2 Liter**  
MOQ: **5 pcs.**

