

moveo®

Medical Workstation

The moveo® Medical Workstation combines ergonomics, flexibility and state-of-the-art technology in a fully integrated solution—designed for versatile use in hospitals, medical practices and care facilities.

With stepless height adjustment and an optional integrated PC, monitor and battery, it offers maximum ease of use for medical applications. Thanks to its modular design, the workstation can be configured individually and adapted to specific requirements.

Whether used as a ward round cart, for creating visit reports, for documentation, monitoring or transporting equipment, the moveo® Medical Workstation adapts perfectly to every task and environment.

Equipment - moveo® Medical Workstation



- 1 support profile, 992 mm, with lockable side door
- 1 shelf, aluminium die-cast, 467 × 349 mm (W × D), with a convenient side release button for stepless height adjustment from 700–1000 mm via gas spring, load capacity up to approx. 18 kg
- 1 shelf, aluminium die-cast, 467 × 349 mm (W × D) and 1 metal drawer, lockable, with ISO modular basket
- Keyboard tray with integrated CHERRY keyboard & mouse pad
- 1 monitor mount, Reell vMX 100, load capacity up to approx. 12 kg
- 1 power strip (3-way)
- 1 multi-purpose clamp
- 4 wall bumpers
- 1 base with 4 electrically conductive twin swivel castors, Ø 125 mm, all 4 lockable
- Cable strain relief for 3 cables
- Central ballast cover for PC
- Isolation transformer box for battery

Dimensions: 576 x 597 x 1364 / 1664 in mm (B x T x H)

Optional: counterweights, isolation transformer, PC, monitor or battery

moveo®

Medical Workstation

Your benefits at a glance

- Ergonomic & flexible: Height-adjustable for comfortable working while seated or standing – suitable for every body height and situation.
- Efficient workflow: Optional integrated IT components and intelligent cable management ensure smooth processes and quick access.
- Long operating time: High-quality batteries provide power for flexible use in different locations.
- Versatile application: Ideal for ward rounds, monitoring, documentation and the safe transport of medical equipment.
- Safe & stable: Robust construction, enclosed surfaces and secure cable routing for maximum stability.
- Durable & future-proof: Premium materials and a modular design for long-term use and easy adaptation.
- Customisable: Available in your corporate colours and with your logo on request – for a professional, consistent appearance.

