

Series 2

Examination- and Reading Lamp




LED

+ highly efficient LED bulb

Colour of lamp	white
Light source	LED, 6,6 W
Light colour	appr. 5.500 K
Light output	appr. 8.000 lux / 500 mm
Light field diameter D10	360 / 720 mm at 500 / 1.000 mm
Life expectancy	appr. 30.000 h (light source)
Electrical mains	220 V, 50 – 60 Hz
Switch	in lamp housing


! Production of lamp in accordance to DIN EN 60601, 3rd Edition



Halogen

+ brilliant, non-reflective illumination of the treatment level

Colour of lamp	white
Light source	halogen bulb with cold light reflector 12 V / 20 W IRC, 10°
Light colour	appr. 3.000 K
Light output	appr. 21.500 lux / 500 mm
Light field diameter D10	140 / 340 mm at 500 / 1.000 mm
Life expectancy	appr. 4.000 h (light source)
Electrical mains	230 V, 50 – 60 Hz
Switch	in lamp housing





Lamp with Double-Joint Articulated Arm

- + Lamp arm with a high quality powder coating and a superior stability against discolouration

Working radius	appr. 1.100 mm
Total weight (net weight)	appr. 2,45 kg
Joint pin	Ø 18, L 35 mm
optional	dimmer on request / international plug / 110 V

with spring-balance technology ¹	without spring-balance technology ²	
L210211A	L210212A	halogen
L210215A	L210216A	LED



Lamp with Flexible Gooseneck Arm

- + Lamp head with an integrated aluminium heat sink for optimal heat dissipation

Working radius	appr. 800 mm
Total weight (net weight)	appr. 2,2 kg
Joint pin	Ø 18, L 35 mm
optional	dimmer on request / international plug / 110 V

L220212A	halogen
L220216A	LED



Reading Lamp

- + night orientation light to find the switch in the dark

Light source	halogen bulb with cold light reflector, 12 V / 20 W, 36°
Life expectancy	appr. 3.000 h (light source)
Switch	in lamp head or cord switch
Working radius	appr. 800 mm
Total weight (net weight)	appr. 1,6 kg / 6,5 kg ¹
Joint pin	Ø 18, L 35 mm
Tube height	appr. 1.070 mm
Roller base	steel, white coated
Base diameter	appr. 635 mm
Castors	twin castors Ø 50 mm, according to DIN EN 60601 electrically conductive and lockable
optional	dimmer on request / international plug / 110 V

L230112A	without roller base
L240112A	without roller base with cord switch
L230122A	on mobile roller base

¹ Stand with Caster Base

Mobile Examination Lamp

- + flexible use: at the right time, at the right place

Working radius	appr. 800 mm
Total weight (net weight)	appr. 7,5 kg
Joint pin	Ø 18, L 35 mm
optional	dimmer on request / international plug / 110 V

LED	halogen	
L220226A	L220222A	fixed height with gooseneck arm
L220236A	L220232A	height adjustable with gooseneck arm
L210226A	L210222A	fixed height with articulated arm
L210236A	L210232A	height adjustable with articulated arm
L210225A	L210221A	fixed height with spring balanced arm
L210235A	L210231A	height adjustable with spring balanced arm

optional: dimmer

optional: handle

