

Branca

made
in
Germany

provita[®]
medical

The ergonomic treatment
unit (patented)

- + Significantly reduced tension on the spine due to an ergonomic seating position
- + Less downtime due to illness / back pain
- + Hygienic work - NEXT TO the bed
- + Foam buffer ensures noiseless handling of the seat mechanism



Branca - now also on **YouTube**

<https://www.youtube.com/watch?v=obP-e3CQxp8>



<https://www.youtube.com/watch?v=CjagW3ghXZU>



provita[®]
medical

provita medical gmbh & co. kg
42929 Wermelskirchen
Germany

Phone: +49 21 93 / 51 05 - 0
Fax: +49 21 93 / 51 05 - 269
info@provita.de
www.provita.de

Branca


The ergonomic treatment
unit (patented)

made
in
Germany

provita[®]
medical



Material:	Stainless steel / Aluminium
Arm- and Leg rest:	Pad of PU-Foam
Saddle seat:	Seat surface of PU-Foam
Castors:	7 twin castors Ø 75 mm 3 out of them with wheel-brake
Seat mechanism:	Aluminium, black powder-coated
Seat height adjustment:	appr. 155 mm with gas spring
Arm-/Leg rest height adjustment:	appr. 830 - 1.100 mm
Horizontal adjusting range:	appr. 190 mm



-  Melamin plate, foldable and removable (appr. 340 × 227 mm)
-  Holder for disposal container and disinfectant bottle up to Ø 69 mm (container and bottle not included in the delivery range)
-  Stainless steel storage basket (appr. 300 × 200 × 100 mm), or optional stainless steel wire-mesh basket (appr. 300 × 200 × 200 mm)