

SHARPS CONTAINER HOLDER - ROUND FOR E.G. "Grathwol"



Sharps Container 0,7 & 1,5 L

Sharps container 0,7 & 1,5 liter
AP Medical DISPO Line

Yellow or red lid

Accessories:

- Wall bracket with T-slot
- Rail clamp with T-slot (DK)
- Rail clamp with T-slot (EU)
- Multibracket with T-slot
- GCX rail bracket with T-slot



Sharps container holder til e.g.
AP Medical 0,7 & 1,5 liter Sharps Containers

Item no.:	JD2550000
Inner Diameter:	Ca. Ø100 mm
Material:	Stainless Steel (AISI 304)
Surface:	Electropolished
Height:	50 mm
Depth:	143 mm (outer dimension)
Thread Diameter:	8 mm
Weight:	189 grams
MOQ:	5 pcs.

SHARPS CONTAINER HOLDER - ROUND FOR E.G. "AP MEDICAL"



Sharps Container 2 Liter

Sharps container 2 litres
AP Medical DISPO Line

Yellow or red lid



Sharps container holder for
1.0-1.3-1.5-2.0-3.0 litre Sharps Containers

Fits perfectly to Klinion EasyCare "Mediq /UD-SON" 1,3 & 2L and 2 liter AP Medical DISPO Line boxes. Yellow or red lid.

Item no.:	JD1490000
Inner Diameter:	Ø 134 mm
Material:	Stainless Steel (AISI 304)
Surface:	Electropolished
Height:	50 mm
Depth:	172 mm (outer dimension)
Thread Diameter:	8 mm
Weight:	220 grams
MOQ:	5 pcs.

SHARPS CONTAINER HOLDER - OVAL FOR E.G. "AP MEDICAL"



Sharps Container 2, 4 & 7 Liter

Sharps container 2 litres
AP Medical DISPO Line

Yellow or red lid

Accessories:
 Wall bracket with T-slot
 Rail clamp with T-slot (DK)
 Rail clamp with T-slot (EU)
 Multibracket with T-slot
 GCX rail bracket with T-slot

Sharps container holder for e.g. AP Medical 2, 4
& 7 litre oval boxes

Fits perfectly to Klinion EasyCare "Mediq /UD-SON" AP Medical DISPO Line boxes.
Yellow or red lid.

Item no.:	JD1590000
Inner Diameter:	133 mm x 133 mm
Material:	Stainless Steel (AISI 304)
Surface:	Electropolished
Height:	50 mm
Depth:	175 mm (outer dimension)
Thread Diameter:	8 mm
Weight:	248 grams
MOQ:	5 pcs.

SHARPS CONTAINER HOLDER - OVAL

For E.G. "AP Medical"



Sharps Container 3 Liter

Sharps container 2 liter AP Medical DISPO Line

Yellow or red lid

Accessories:

Wall bracket with T-slot
 Rail clamp with T-slot (DK)
 Rail clamp with T-slot (EU)
 Multibracket with T-slot
 GCX rail bracket with T-slot

Sharps container holder for e.g. AP Medical
 3 liter, oval container

Item no.: **JD1690000**
 Inner Diameter: **143 mm x 143 mm**
 Material: **Stainless Steel (AISI 304)**
 Surface: **Electropolished**
 Height: **50 mm**
 Depth: **175 mm (outer dimension)**
 Thread Diameter: **8 mm**
 Weight: **319 grams**
 MOQ: **5 pcs.**

SHARPS CONTAINER HOLDER - SQUARED

For E.G. "Sharpsafe"



Sharps Container 4 & 7 Liter

Sharps containers 4 / 7 liter
 "Sharpsafe" container

Yellow or red lid

Sharps container holder for e.g. "Sharpsafe"
 4 & 7 liter squared container

Item no.: **JD1540000**
 Inner Diameter: **156 mm x 156 mm**
 Material: **Stainless Steel (AISI 304)**
 Surface: **Electropolished**
 Height: **50 mm**
 Depth: **208 mm (outer dimension)**
 Thread Diameter: **8 mm**
 Weight: **280 grams**
 MOQ: **5 pcs.**



SHARPS CONTAINER HOLDER - SQUARED

For E.G. "Sharpsafe"



Sharps Container 7,5 Liter

Sharps container "Sharpsafe"
7.5 liter QUIVER

Yellow or red lid

Sharps container holder for e.g.
"Sharpsafe Quiver"
7,5 liter, squared container (long model)

Item no.: **JD2790000**
Inner Diameter: **101 mm x 147 mm**
Material: **Stainless Steel (AISI 304)**
Surface: **Electropolished**
Height: **50 mm**
Depth: **153 mm (outer dimension)**
Thread Diameter: **8 mm**
Weight: **247 grams**
MOQ: **5 pcs.**

Accessories:
Wall bracket with T-slot
Rail clamp with T-slot (DK)
Rail clamp with T-slot (EU)
Multibracket with T-slot
GCX rail bracket with T-spur

SHARPS CONTAINER HOLDER - SQUARED

For E.G. "Europlásticos"



Sharps Container 1 Liter

Sharps container for 1 liter boks
from Europlásticos

Yellow or red lid

Sharps container holder for e.g. "Sharpsafe"
1 liter squared container

Item no.: **JD1700000**
Inner Diameter: **96 mm x 96 mm**
Material: **Stainless Steel (AISI 304)**
Surface: **Electropolished**
Height: **50 mm**
Depth: **208 mm (outer dimension)**
Thread Diameter: **8 mm**
Weight: **200 grams**
MOQ: **5 pcs.**



SHARPS CONTAINER HOLDER - SQUARED

Other Assortment - e.g. "Safebox"



Sharps Container 4L, 6L & 11 Liter

Sharps Container "Safebox" 4L, 6L & 11L

Yellow or red lid

Accessories:

Wall bracket with T-slot
 Rail clamp with T-slot (DK)
 Rail clamp with T-slot (EU)
 Multibracket with T-slot
 GCX rail bracket with T-slot



Sharps Container Holder for e.g.
 "Sharpsafe Quiver"
 7,5 liter, squared container (long model)

Item no.: **JD2790000**
 Inner Diameter: **101 mm x 147 mm**
 Material: **Stainless Steel (AISI 304)**
 Surface: **Electropolished**
 Height: **50 mm**
 Depth: **153 mm (outer dimension)**
 Thread Diameter: **8 mm**
 Weight: **247 grams**
 MOQ: **5 pcs.**

SHARPS CONTAINER HOLDER - ROUND

Other Assortment



Sharps Container Holder for
 0,5- & 1 liter round container

Fits e.g. Dispo 0,5 litre box from Grathwol.

Item no.: **JD1500000**
 Inner Diameter: **Ø73 mm**
 Material: **Stainless Steel (AISI 304)**
 Surface: **Electropolished**
 Height: **50 mm**
 Depth: **112 mm (outer dimension)**
 Thread Diameter: **8 mm**
 Weight: **150 grams**
 MOQ: **5 pcs.**



Sharps Container Holder for
 0.6 litre - round container

Fits e.g. Multisafe 1 l. & USON 500 ml.

Item no.: **JD1510000**
 Inner Diameter: **Ø89 mm**
 Material: **Stainless Steel (AISI 304)**
 Surface: **Electropolished**
 Height: **50 mm**
 Depth: **134 mm (outer dimension)**
 Thread Diameter: **8 mm**
 Weight: **171 grams**
 MOQ: **5 pcs.**

SHARPS CONTAINER HOLDER - ROUND

Other Assortment



**Sharps Container Holder for e.g. Mediq /USON
5 & 7 litre round container**

Item no.: **JD2640000**
 Inner Diameter: **Ø200 mm**
 Material: **Stainless Steel (AISI 304)**
 Surface: **Electropolished**
 Height: **50 mm**
 Depth: **244 mm (outer dimension)**
 Thread Diameter: **8 mm**
 Weight: **286 grams**
 MOQ: **5 pcs.**



**Sharps Container Holder for e.g. USON 11 litres
Jokey container from Mediq**

Item no.: **JD2660000**
 Inner Diameter: **Ø251 mm**
 Material: **Stainless Steel (AISI 304)**
 Surface: **Electropolished**
 Height: **50 mm**
 Depth: **295 mm (outer dimension)**
 Thread Diameter: **8 mm**
 Weight: **345 grams**
 MOQ: **5 pcs.**



**Sharps Container Holder for e.g. USON 21 litre
and Jokey container from Mediq**


Item no.: **JD2670000**
 Inner Diameter: **Ø280 mm**
 Material: **Stainless Steel (AISI 304)**
 Surface: **Electropolished**
 Height: **50 mm**
 Depth: **324 mm (outer dimension)**
 Thread Diameter: **8 mm**
 Weight: **368 grams**
 MOQ: **5 pcs.**



Accessories:
 Wall bracket with T-slot
 Rail clamp with T-slot (DK)
 Rail clamp with T-slot (EU)
 Multibracket with T-slot
 GCX rail bracket with T-slot

WIRE BASKETS & WIRE PRODUCTS

Stainless steel



**Wire basket 6 litres, conical
with hook for T-slot**

**Wire basket 7 litres, conical
with hook for T-slot**

**Wire basket 15 litres, conical
with hook for T-slot**

Item no.: **JD1610000**
 Dimensions: **Width x Depth x Height
231mm x 161mm x 242mm
Stainless Steel (AISI 304)**
 Material: **Electropolished**
 Surface:
 Weight: **853 grams**
 Volumen: **6 Liter**
 MOQ: **5 pcs.**

Accessories for all 3 types of wire baskets:

Wall bracket with T-slot
 Rail clamp with T-slot (DK)
 Rail clamp with T-slot (EU)
 Multibracket with T-slot
 GCX rail bracket with T-spur

Item no.: **JD1610101**
 Dimensions: **Width x Depth x Height
623mm x 217mm x 132mm
Stainless Steel (AISI 304)**
 Material: **Electropolished**
 Surface:
 Weight: **667 grams**
 Volumen: **7 Liter**
 MOQ: **5 pcs.**

Item no.: **JD1610201**
 Dimensions: **Width x Depth x Height
337mm x 233mm x 267mm
Stainless Steel (AISI 304)**
 Material: **Electropolished**
 Surface:
 Weight: **959 grams**
 Volumen: **15 Liter**
 MOQ: **5 pcs.**